

# Omni-C™ Technology


*Delivers SNPs to Genome Conformation In a Single Assay*

Powered by a sequence-independent endonuclease, Omni-C™ technology represents the next-evolution in Hi-C assays delivering unbiased long-range reads and highly uniform sequence coverage.

Gain a comprehensive view of the genome from SNPs to indels, SNVs to phasing, and loops to topologically associated domains (TADs), all from a single NGS library.

## ACHIEVE SHOTGUN-LIKE COVERAGE

Omni-C libraries contain uniform, shotgun-like coverage enabling downstream applications sensitive to variant detection.



Coverage distribution of different Hi-C libraries (colored lines) compared to a shotgun library (dotted line).

## CAPTURE LONG-RANGE INFORMATION

Omni-C libraries offer greater long-range information compared to other Hi-C approaches. The result is greater support for chromatin 3-D interactions such as TADs and loops enabling comprehensive genome conformation mapping.

**A.**


| Library Type   | % cis <1 kbp | % cis >1 kbp | % unique molecules at 300 M read pairs |
|----------------|--------------|--------------|--|
| Omni-C         | 5%           | 95%          | 73%                                    |
| Hi-C Multi-RE  | 20%          | 80%          | 69%                                    |
| Hi-C Single-RE | 34%          | 66%          | 50%                                    |

Comparison of Omni-C libraries with multi-restriction enzyme (RE) and single-RE Hi-C libraries.

A. Proportion of long-range reads for each library type

B. Contact matrices at 4 kb resolution from 800 M total reads per library. Loops are circled with the number of raw reads supporting each contact indicated.

**B.**



To place an order or for more information:

visit us at [www.dovetailgenomics.com](http://www.dovetailgenomics.com) or send an email to [info@dovetail-genomics.com](mailto:info@dovetail-genomics.com)

## GENETIC VARIANTS AND PHASED GENOMES FROM A SINGLE LIBRARY

The shotgun-like characteristics of Omni-C data enable genotyping and haplotype phasing.

| Library Type  | True Positive    | False Positive | False Negative | Sensitivity  | Precision    |
|---------------|------------------|----------------|----------------|--------------|--------------|
| Shotgun       | 2,696,291        | 9,270          | 6,814          | 99.7%        | 99.7%        |
| <b>Omni-C</b> | <b>2,666,339</b> | <b>20,081</b>  | <b>36,766</b>  | <b>98.6%</b> | <b>99.3%</b> |
| RE Based Hi-C | 2,387,235        | 33,554         | 315,870        | 88.3%        | 98.6%        |


*Homozygous SNP calls sampled at 300 M read pairs.*

| Library Type  | #Variants Phased | % Heterozygous SNPs phased | Switch Error Rate | Largest Phase Block Size | #Chromosomes phased end to end |
|---------------|------------------|----------------------------|-------------------|--------------------------|--------------------------------|
| Shotgun       | 2,229,492        | 81.44%                     | 0.0036            | 28.2 kb                  | 0                              |
| <b>Omni-C</b> | <b>2,299,248</b> | <b>84.00%</b>              | <b>0.0100</b>     | <b>248.01 Mb</b>         | <b>23</b>                      |
| RE Based Hi-C | 1,986,467        | 72.60%                     | 0.0357            | 247.95 Mb                | 23                             |

*GM12878 haplotype assembly metrics using high confidence heterozygous SNPs from the Illumina Platinum Genome.*

## EXPLORE THE LANDSCAPE OF STRUCTURAL VARIANTS

Omni-C libraries capture large structural (>1 Mbp) variant information such as indels, inversions and chromosomal translocations.




*Validated large SVs found in breast cancer cell line HCC1187 representing an inversion, a deletion and inter-chromosomal translocation*

## EXTRACT MORE DATA FROM EVERY SEQUENCING RUN

Omni-C technology captures genome-wide topology at single nucleotide resolution in a single library prep.

Therefore, in addition to being ideal for the study of 3-D genome conformation, it is suitable for applications traditionally addressed by whole genome sequencing.



To place an order or for more information:

visit us at [www.dovetailgenomics.com](http://www.dovetailgenomics.com) or send an email to [info@dovetail-genomics.com](mailto:info@dovetail-genomics.com)